

#### **Missouri Artificial Vagina & Accessories**

• Kit includes: 16" or 22" latex A/V, leather jacket, two plastic 250 cc collection bottles, bottle adapter, semen filters, filling system, A/V conditioning solution and ten disposable inner liners

HAR-VET™

- Latex A/V Liner, 16" or 22"
- A/V filling valve air and water
- One-piece latex liner with tapered cone
- Leather jacket and reinforced handle assist in collection
- Fills with 3/8" Filling System or the E-Z A/V Filler
- Regular use of Conditioning Solution for A/V Liners will reduce bladder sticking and help prevent the growth of mold and mildew
- Approximate weight when full is 12 lbs

| Item  | Product Description  | UOM    |
|-------|--|--------|
| 61629 | 22" Missouri A/V Complete Kit  | Each   |
| 61622 | Leather Jacket for 22" Missouri A/V                                  | Each   |
| 61624 | Latex Liner for 22" Missouri A/V                                     | Each   |
| 17145 | Disposable A/V Liner for 22" Missouri A/V                            | 10/Bag |
| 17104 | Insulated Neoprene A/V Cover for 22" Missouri A/V                    | Each   |
| 61630 | 16" Missouri A/V Complete Kit  | Each   |
| 61625 | Latex Liner 16" Missouri A/V   | Each   |
| 61623 | Leather Jacket for 16" Missouri A/V                                  | Each   |
| 61636 | Disposable A/V Liner for 16" Missouri A/V                            | 10/Bag |
| 17105 | Insulated Neoprene A/V Cover for 16" Missouri A/V                    | Each   |
| 17120 | Insulated Bottle Cover for Missouri A/V                              | Each   |
| 17135 | 3/8" Filling System for Missouri A/V                                 | Each   |
| 61626 | Semen Collection Bottle Missouri / Har-Vet <sup>™</sup> A/V - 250 cc | 2/Bag  |
| 61631 | Bottle Adapter for Missouri A/V                                      | Each   |
| 61642 | A/V Semen Filter - 6" x 2-1/4"                                       | 25/Bag |
| 17140 | E-Z A/V Filler   | Each   |

#### Storage & Care of Missouri A/V Liners

- NEVER expose to petroleum-based products
- Do not stretch or fold during storage
- Store in a COOL, DARK, DRY location

#### **Cleaning of Latex Rubber**

- DO NOT use harsh detergents, antibacterial or degreaser soaps. Use non-ionic soap
- RINSE WELL
- DRY latex rubber immediately! Remove all water from the bladder. Hang the liner to dry. Remove water stem while drying

Equine A/V is guaranteed free from defects for 90 days from the ship date.

Reproductior

#### Colorado Artificial Vagina & Accessories

- Kit includes: A/V housing, water jacket liner, ten 7" strong securing bands, four plastic 250 cc collection bottles, 25 semen filters, water filling system and ten disposable inner liners
- 22-1/2" rigid one-piece plastic housing with handle
- Housing has a large 6" inner diameter
- A/V filling valve air and water


•

- Fills with the 1/2" Filling System or the E-Z A/V Filler
- Approximate weight when full is 25 lbs

| Item  | Product Description  | UOM    |
|-------|--|--------|
| 61645 | Colorado A/V Complete Kit  | Each   |
| 61644 | Housing Only for Colorado A/V  | Each   |
| 61646 | Latex Water Jacket Liner for Colorado A/V                            | Each   |
| 17035 | Latex Collection Liner with Flared End for Colorado A/V              | Each   |
| 17137 | 1/2" Filling System for Colorado A/V                                 | Each   |
| 61648 | 7" Rubber Band for Colorado A/V                                      | 10/Bag |
| 61649 | Disposable A/V Liner without Filter for Colorado A/V                 | 10/Bag |
| 61642 | A/V Semen Filter - 6" x 2-1/4"                                       | 25/Bag |
| 17140 | E-Z A/V Filler   | Each   |
| 17121 | Protective Cone - Colorado A/V                                       | Each   |
| 61626 | Semen Collection Bottle Missouri / Har-Vet <sup>™</sup> A/V - 250 cc | 2/Bag  |









#### Har-Vet<sup>™</sup> Artificial Vagina & Accessories

- Kit includes: A/V housing, one latex water jacket liner, two sets of strong securing bands, two plastic 250 cc collection bottles, 25 semen filters, ten disposable inner liners and instruction manual
- Similar size to the Japanese A/V; housing inside diameter 4-3/4"
- 22-1/2" rigid one-piece plastic housing with handle
- Metal filling port allows for rapid filling
- Approximate weight when full is 8 lbs

| ltem  | Product Description  | UOM    |
|-------|--|--------|
| 61602 | Har-Vet <sup>™</sup> A/V Complete Kit                                | Each   |
| 61600 | Har-Vet <sup>™</sup> A/V Housing Only                                | Each   |
| 61603 | Har-Vet <sup>™</sup> Latex Water Jacket Liner                        | Each   |
| 61606 | Har-Vet <sup>™</sup> Disposable A/V Liner - Open End                 | 10/Bag |
| 61609 | Har-Vet <sup>™</sup> Rubber Band - Large                             | 2/Bag  |
| 61610 | Har-Vet <sup>™</sup> Rubber Band - Small                             | 2/Bag  |
| 61626 | Semen Collection Bottle Missouri / Har-Vet <sup>™</sup> A/V - 250 cc | 2/Bag  |
| 61642 | A/V Semen Filter - 6" x 2-1/4"                                       | 25/Bag |
| 61607 | Whirl Pak <sup>®</sup> Collection Bag, 18 oz                         | 10/Bag |



#### Disposable Artificial Vagina Kit

- Kit includes: 20" outer housing, inner vinyl liner, four filters, one plastic 250 cc collection bottle, four cable ties and four clamps to close water filling stem
- Housing inside diameter 6"
- Inner bladder is made of soft vinyl and is attached to the housing
- Inner bladder is filled with water
- Ideal for limited use in clinics or stables
- Collection same as a Colorado, Har-Vet<sup>™</sup>, Missouri or Japanese A/V
- Designed for single use to prevent the possibility of disease transfer
- Single use item

| Item  | Product Description       | UOM  |
|-------|---------------------------|------|
| 61640 | Disposable Equine A/V Kit | Each |

#### Pony / Mini Artificial Vagina Complete Kit

- Kit includes: rubber A/V housing, leather jacket with handle, one latex water liner, ten disposable liners with filters, two strong securing bands and a water filling system
- A/V filling valve air and water

**Global Locations:** 

- Fills with the 3/8" Filling System or the E-Z A/V Filler
- Housing is 12" long with a 4" inside diameter
- For use with ponies and miniature stallions

| Item  | Product Description  | UOM    |
|-------|--|--------|
| 17026 | 12" Pony / Mini A/V Complete Kit                                     | Each   |
| 17106 | Leather Jacket for Pony / Mini A/V                                   | Each   |
| 17107 | Latex Inner Liner for Pony / Mini A/V                                | Each   |
| 17109 | Disposable A/V Liner with Filter for Pony / Mini A/V                 | 10/Bag |
| 17111 | Disposable A/V Liner without Filter for Pony / Mini A/V              | 10/Bag |
| 61626 | Semen Collection Bottle Missouri / Har-Vet <sup>™</sup> A/V - 250 cc | 2/Bag  |
| 61642 | A/V Semen Filter - 6" x 2-1/4"                                       | 25/Bag |
| 17135 | 3/8" Filling System for Missouri A/V                                 | Each   |
| 17140 | E-Z A/V Filler   | Each   |







#### Standard Bovine Artificial Vagina Complete

- Complete A/V includes: optional size housing fitted with a screw cap water filling valve, one latex inner liner, two strong securing bands, one latex director cone and one graduated semen collection vial
- Housing made of high quality black rubber tubing with 2-3/4" inside diameter

| Item  | Product Description   | UOM  |
|-------|---|------|
| 14706 | Standard 14-1/2" Bovine A/V Complete                        | Each |
| 14701 | Standard 14-1/2" Bovine A/V Housing Only with Filling Valve | Each |
| 14705 | Standard 16-1/2" Bovine A/V Complete                        | Each |
| 14700 | Standard 16-1/2" Bovine A/V Housing Only with Filling Valve | Each |



- Complete A/V includes: 16-1/2" housing, one latex inner liner, two strong securing bands, one latex director cone, one graduated semen collection vial and one water filling valve
- Housing made of molded rubber with optimal rigidity
- Smooth finish inside and out
- Rounded ends allow easy cleaning
- 5/8" water filling hole fitted with screw cap water filling valve

| Item  | Product Description                                       | UOM  |
|-------|---|------|
| 14709 | Deluxe 16-1/2" Bovine A/V Complete                        | Each |
| 14708 | Deluxe 16-1/2" Bovine A/V Housing Only with Filling Valve | Each |
| 14707 | Deluxe 16-1/2" Bovine A/V Housing Only                    | Each |

#### Bovine Artificial Vagina Accessories (MGI)

| Item  | Product Description                               | UOM     |
|-------|---|---------|
| 14610 | Inner Liner, Smooth Surface Latex - 23"           | Each    |
| 14651 | Semen Collection Cone, Latex - 10" without Vent   | Each    |
| 14655 | Semen Collection Cone, Disposable Plastic - 11"   | 100/Box |
| 50014 | Semen Collection Cone, Disposable Plastic - 14"   | 25/Bag  |
| 14657 | Disposable A/V Liner, 20" Bovine                  | 10/Bag  |
| 17140 | E-Z A/V Filler                                    | Each    |
| 17142 | Economy A/V Filler                                | Each    |
| 49101 | A/V Disposable Collection Centrifuge Vial - 15 ml | 100/Box |
| 14800 | A/V Filling Valve - Air & Water                   | Each    |
| 14710 | A/V Filling Valve Screw Cap                       | Each    |
| 15800 | A/V Air Filling Jug Pump                          | Each    |



#### Ovine Artificial Vagina Complete & Accessories

 Complete A/V includes: 8" rubber A/V housing with a 2-1/4" inside diameter, air and water valve, one latex inner liner, one latex collection liner, two strong securing bands, one graduated semen collection vial and one jug pump

| Item  | Product Description                       | UOM  |
|-------|---|------|
| 15010 | Sheep A/V Complete                        | Each |
| 15005 | Sheep A/V Housing Only with Filling Valve | Each |
| 15004 | Sheep A/V Housing Only                    | Each |
| 15011 | Sheep A/V Inner Liner                     | Each |



www.maianimalhealth.com









#### "The Tickler"

- Includes Electronic Control Unit fitted with vinyl protective jacket, Standard Bull Probe, Ram Probe, SLA battery charger and carrying case
- 12-Volt battery lead available
- Various probes available bull or ram

| ltem  | Product Description  | UOM  |
|-------|--|------|
| 14900 | "The Tickler" <sup>™</sup> - Electronic Control Unit (EJ6CR7) - CGS Ejaculator Kit | Each |
| 14901 | 12V Battery Lead   | Each |
| 14902 | 12V US Battery Charger   | Each |
| 14903 | Ram Probe  | Each |
| 14907 | Standard Bull Probe  | Each |
| 14905 | Dual Mode Bull Probe   | Each |
| 14909 | Multi Mode Bull Probe  | Each |





| ltem  | Product Description              | UOM  |
|-------|----------------------------------|------|
| 17100 | A/V Conditioning Solution - 4 oz | Each |

## Artificial Vagina Semen Filter

- Fits any A/V using a baby bottle for collection including Har-Vet<sup>™</sup>, Missouri, Colorado, • Disposable and Pony A/V
- Used to remove debris, non-semen or seminal fluids from the ejaculate, leaving whole semen for laboratory evaluation and A.I. procedures
- Low cost
- Single use item

| ltem  | Product Description            | UOM    |
|-------|--------------------------------|--------|
| 61642 | A/V Semen Filter - 6" x 2-1/4" | 25/Bag |









•

.



| Semen ( | Sollection Bottle  |       |
|---------|--|-------|
| ltem    | Product Description  | UOM   |
| 61626   | Semen Collection Bottle Missouri / Har-Vet <sup>™</sup> A/V - 250 cc | 2/Bag |







#### E-Z Artificial Vagina Filler

- Fits onto hose attachment
- Easy to use without getting wet
- Quick release connector for easy attachment
- Brass fittings
- Fits most A/Vs

| Item  | Product Description | UOM  |
|-------|---------------------|------|
| 17140 | E-Z A/V Filler      | Each |

#### Economy Artificial Vagina Filler

| ltem  | Product Description | UOM  |
|-------|---------------------|------|
| 17142 | Economy A/V Filler  | Each |

## Hanger for Latex Liner

| ٠ | For use in properly drying and storing latex liners |
|---|---|
|   |   |

| ltem  | Product Description    | UOM  |
|-------|------------------------|------|
| 17101 | Hanger for Latex Liner | Each |

## Valve for Air / Liquid Core Housing

• Replacement filling valve for Missouri, Colorado and Mini A/V

| Item  | Product Description                 | UOM  |
|-------|-------------------------------------|------|
| 17800 | Valve for Air / Liquid Core Housing | Each |

#### Stallion Condom

- Latex
- 12-3/4" long, 5-3/4" wide at base and 3-3/8" wide at top opening
- .0125" wall thickness at open end

| Item  | Product Description | UOM  |
|-------|---------------------|------|
| 61650 | Stallion Condom     | Each |

## Alconox<sup>®</sup> Detergent Powder

- Concentrated detergent for removing organic contaminants
- May be used on soft metals and latex rubber

| Item  | Product Description                          | UOM  |
|-------|--|------|
| 17103 | Alconox <sup>®</sup> Detergent Powder, 4 lbs | Each |









**Artificial Vaginas & Accessories** 







Europe: +44 (0) 1225 791 184

Canine - Semen Collection & Tie-Up<sup>™</sup>

#### Canine Semen Collection & Artificial Insemination Kit (MGI)

Complete kit for the collection of canine semen. Kit includes one each:

- Latex collection cone and 15 cc graduated vial •
- Disposable 5 cc syringe
- 4" insemination pipette for small dogs
- 9" insemination pipette for large dogs

| - 0   | insemination pipette for large dogs                   |      |
|-------|---|------|
| ltem  | Product Description                                   | UOM  |
| 19250 | Canine Semen Collection & Artificial Insemination Kit | Each |
|       |   |      |

Pair of disposable gloves

Single use item

Package of lubrication jelly

#### (MQi) Canine Disposable A/V Liner Kit & Collection Cones

A/V Liner Kit includes a disposable liner made of 2.5 mil polyethylene film attached to a 15 cc graduated vial. It may be used as a collection cone or as a liner inside a collection cone.

• Single use item

|       | •  |         |
|-------|--|---------|
| ltem  | Product Description                                      | UOM     |
| 50015 | Canine Disposable A/V Liner Kit                          | 10/Bag  |
| 50014 | Semen Collection Cone, Disposable Plastic - 14", 2.5 mil | 25/Bag  |
| 14655 | Semen Collection Cone, Disposable Plastic - 11", 3 mil   | 100/Box |
| 14650 | Semen Collection Cone, Latex - 10" with Vent             | Each    |
| 14651 | Semen Collection Cone, Latex - 10" without Vent          | Each    |




#### (MQi) PATENTED Tie-Up<sup>™</sup> Canine Intravaginal A.I. Kits

- Soft, plastic coated sponge simulates male penis, triggering normal physiological events • in the female canine that would occur during actual breeding
- After insertion, female's vaginal contractions close down on sponge simulating actual "tying up"
- A.I. catheter runs through center of vaginal delivery device. Proximal end has non-traumatic STIP<sup>™</sup> Tip and distal end has flex tip for syringe or cap attachment
- Vaginal contractions on sponge hold sponge in place and stop semen leakage; the contractions aspirate the semen into the uterine body cavity simulating the natural act
- Kit includes: canine disposable A/V liner, 15 cc graduated vial and Canine Intravaginal A.I. Device
- Three sizes available: 2", 4" and 6" individually or in kits
- Available non-sterile or gamma irradiated
- Single use item .

**Global Locations:** 

| Item    | Product Description  | UOM   |
|---------|--|-------|
| 18200-2 | Tie-Up <sup>™</sup> Canine Intravaginal A.I. Kit, 2"                       | Each  |
| 18200-4 | Tie-Up <sup>™</sup> Canine Intravaginal A.I. Kit, 4"                       | Each  |
| 18200-6 | Tie-Up <sup>™</sup> Canine Intravaginal A.I. Kit, 6"                       | Each  |
| 18206-2 | Tie-Up <sup>™</sup> Canine Intravaginal A.I. Device, 2" - Gamma Irradiated | 5/Bag |
| 18206-4 | Tie-Up <sup>™</sup> Canine Intravaginal A.I. Device, 4" - Gamma Irradiated | 5/Bag |
| 18206-6 | Tie-Up <sup>™</sup> Canine Intravaginal A.I. Device, 6" - Gamma Irradiated | 5/Bag |
| 18205-2 | Tie-Up <sup>™</sup> Canine Intravaginal A.I. Device, 2" - Non-Sterile      | 5/Bag |
| 18205-4 | Tie-Up <sup>™</sup> Canine Intravaginal A.I. Device, 4" - Non-Sterile      | 5/Bag |
| 18205-6 | Tie-Up <sup>™</sup> Canine Intravaginal A.I. Device, 6" - Non-Sterile      | 5/Bag |







#### Infusion Pipettes

Single use pipette, individually packaged for canine artificial insemination. Pipette is 4" or 9" long

- with flex tip adapter to fit a syringe.
  - Individually packaged
  - Gamma irradiated or EtO exposed
  - Single use item

| ltem  | Product Description                                   | UOM    |
|-------|---|--------|
| 26600 | Infusion Pipette, 4" with Flex Tip - Gamma Irradiated | 12/Bag |
| 19020 | Infusion Pipette, 9" with Flex Tip - Gamma Irradiated | 25/Bag |



Single use pipette with STIP<sup>™</sup> Tip smooth bulb end allows for less traumatic passage of pipette without complication. Tapered tip prevents back vacuum of semen when withdrawing the pipette. Narrow outlet port is attached to a 4" or 10" pipette to deliver semen further into anterior vagina. Flex tip adapter fits onto syringe.

- Individually packaged
- Gamma irradiated or EtO exposed
- Single use item

| Item  | Product Description                                    | UOM    |
|-------|--|--------|
| 50405 | STIP <sup>™</sup> Tip Canine Insemination Pipette, 4"  | 12/Bag |
| 50400 | STIP <sup>™</sup> Tip Canine Insemination Pipette, 10" | 25/Bag |

#### Canine Bullet Tip Culture Swab

- Designed specifically for use in canines
- 14" double-guarded swab consists of an outer soft PVC sleeve with 'X' perforated tip covering an inner tubular guard which houses the calcium alginate swab style collection tip
- Assists in preventing outside contaminants from affecting the sample
- Individually packaged
- Gamma irradiated or EtO exposed
- Single use item

| Item  | Product Description            | UOM    |
|-------|--------------------------------|--------|
| 50020 | Canine Bullet Tip Culture Swab | 25/Bag |



Kenny Skim Milk Canine Semen Extenders provide the user with the following advantages:

- Prevents agglutination of spermatozoa so forward motility can be determined
- Canine semen extender with gentamicin helps destroy infectious pathogens
- Buffers the pH of antibiotics added to the semen
- Maintains a higher percent of forward motility over time as compared to raw semen
- Increases the insemination dose volume making it easier to handle and minimizing
- spermatozoa loss in equipment such as insemination pipettes
- Allows for storage of cooled semen for use up to 48 hours after collection

| Item  | Product Description  | UOM  |
|-------|--|------|
| 50420 | Canine Semen Extender without Antibiotics - Non-Returnable | Each |
| 50423 | Canine Semen Extender with Gentamicin - Non-Returnable     | Each |











#### **INRA 96 Semen Extender**

- Pre-mixed ready to use liquid
- . Ideal for problem stallions
- Contains Gentamycin, Penicillin and Amphotericin B
- Store refrigerated for up to one week after opening
- Sterile

| ltem  | Product Description                             | UOM  |
|-------|---|------|
| 17553 | INRA 96 Semen Extender, 200 ml - Non-Returnable | Each |



#### Har-Vet<sup>™</sup> Semen Extenders & Accessories • Prevents agglutination of spermatozoa so forward motility can be determined • Helps destroy infectious pathogens when an antibiotic has been added • Buffers the pH of antibiotics added to the semen

- Maintains a higher percent of forward motility over time as compared to raw semen
- Increases the insemination dose volume making it easier to handle and minimizing spermatozoa loss in equipment such as insemination pipettes
- Allows for storage of cooled semen for use up to 48 hours after collection
- Semen extender and sterile-filtered water available as a kit or sold separately
- Packaged in moistureproof foil to ensure stable and long lasting shelf life •

| Item  | Product Description   | UOM    |
|-------|---|--------|
| 01500 | Har-Vet <sup>™</sup> Semen Extender Kit with Gentamicin Sulfate - Non-Returnable                        | Each   |
| 01501 | Har-Vet <sup>™</sup> Semen Extender with Gentamicin Sulfate - Powder Only - Non-Returnable              | Each   |
| 01502 | Har-Vet <sup>™</sup> Semen Extender Kit with Amikacin & Potassium Penicillin - Non-Returnable           | Each   |
| 01503 | Har-Vet <sup>™</sup> Semen Extender with Amikacin & Potassium Penicillin - Powder Only - Non-Returnable | Each   |
| 01510 | Har-Vet <sup>™</sup> Semen Extender Kit without Antibiotics - Non-Returnable                            | Each   |
| 01511 | Har-Vet <sup>™</sup> Semen Extender without Antibiotics - Powder Only - Non-Returnable                  | Each   |
| 01520 | Har-Vet <sup>™</sup> Semen Extender Kit with Ticarcillin - Non-Returnable                               | Each   |
| 01521 | Har-Vet <sup>™</sup> Semen Extender with Ticarcillin - Powder Only - Non-Returnable                     | Each   |
| 01525 | 100 ml Bottled Water for Semen Extender Powder - Sterile-Filtered                                       | 25/Box |
| 01526 | 92 ml Bottled Water for Semen Extender Powder - Sterile-Filtered  | 25/Box |

S. M

HAR-VET™



8





#### DVM Rapid Test<sup>™</sup> I

| Product Description  | UOM   |
|--|---|
| Sperm Concentration Test - Bovine, Equine & Porcine - for Original DVM Rapid Test <sup>™</sup> I | 25/Box  |
| Sperm Concentration Test - Bovine, Equine & Porcine - for Original DVM Rapid Test <sup>®</sup> I | 125/Box   |
| Sperm Concentration Test - Bovine, Equine & Porcine - for Original DVM Rapid Test <sup>™</sup> I | 1,000/Box   |
| Sperm Concentration Reagent - 25 ml - for Original DVM Rapid Test <sup>™</sup> I                 | Each  |
| Sperm Concentration Reagent - 100 ml - for Original DVM Rapid Test <sup>™</sup> I                | Each  |
| Primary Turbidity Standards Kit - for Original DVM Rapid Test <sup>™</sup> I                     | Each  |
| Sperm Concentration Test - Canine - for Original DVM Rapid Test <sup>™</sup> I                   | 25/Box  |
|  | Sperm Concentration Test - Bovine, Equine & Porcine - for Original DVM Rapid Test <sup>™</sup> I<br>Sperm Concentration Test - Bovine, Equine & Porcine - for Original DVM Rapid Test <sup>™</sup> I<br>Sperm Concentration Test - Bovine, Equine & Porcine - for Original DVM Rapid Test <sup>™</sup> I<br>Sperm Concentration Reagent - 25 ml - for Original DVM Rapid Test <sup>™</sup> I<br>Sperm Concentration Reagent - 100 ml - for Original DVM Rapid Test <sup>™</sup> I<br>Primary Turbidity Standards Kit - for Original DVM Rapid Test <sup>™</sup> I |



- Portable, turn key and battery-operated handheld analyzer kit
- Lightweight portable case contains all components required to perform test
- Low cost and simple to operate
- All tests may be performed in 10 minutes or less
- Proven accuracy
- Simple, one-step test procedure
- Direct reading of results; no charts needed
- Test kits sold separately
- Optional infrared router for transfering data to computer
- DVM Rapid Test<sup>™</sup> II Multi-Test Analyzer Tests:
  - Serum / Plasma IgG Bovine, Equine and Camelid
  - Colostrum IgG Bovine and Equine
  - Sperm Concentration Equine
- DVM Rapid Test<sup>™</sup> II Sperm Concentration Analyzer Tests:
  - Sperm Concentration Equine, Bovine, Porcine and Canine

| Item    | Product Description  | UOM  |
|---------|--|------|
| 49030   | DVM Rapid Test <sup>™</sup> II - Multi-Test Analyzer                 | Each |
| 49035   | DVM Rapid Test <sup>™</sup> II - Sperm Concentration Analyzer        | Each |
| 49586-1 | Primary Turbidity Standards Kit - For DVM Rapid Test <sup>™</sup> II | Each |

#### IgG Serum / Plasma & Colostrum Test Kits - Equine

| Item    | Product Description  | UOM    |
|---------|--|--------|
| 49500-1 | IgG Serum / Plasma Test Kit - Equine - Non-Returnable, Refrigerate | 12/Box |
| 49505-1 | IgG Colostrum Test Kit - Equine - Non-Returnable, Refrigerate      | 12/Box |

#### IgG Serum / Plasma & Colostrum Test Kits - Bovine

| Item    | Product Description  | UOM    |
|---------|--|--------|
| 49550-1 | IgG Serum / Plasma Test Kit - Bovine - Non-Returnable, Refrigerate | 12/Box |
| 49555-1 | IgG Colostrum Test Kit - Bovine - Non-Returnable, Refrigerate      | 12/Box |

#### IgG Serum / Plasma Test Kit - Camelid

| Item    | Product Description   | UOM    |
|---------|---|--------|
| 49530-1 | IgG Serum / Plasma Test Kit - Camelid - Non-Returnable, Refrigerate | 12/Box |

#### **Sperm Concentration Test - Equine**

| Item    | Product Description   | UOM    |
|---------|---|--------|
| 49510-1 | Sperm Concentration Test - Equine, Bovine, Porcine & Canine | 25/Box |

#### DVM Rapid Test<sup>™</sup> II Router

| Item  | Product Description                   | UOM  |
|-------|---------------------------------------|------|
| 49034 | DVM Rapid Test <sup>™</sup> II Router | Each |



\* Test Kit and Pipettes sold separately.

#### Tip:

Use the DVM Rapid Test<sup>™</sup> II Router to easily transfer data from the DVM Rapid Test<sup>™</sup> II to a computer!









Mini Item

# Pipettors & Accessories

Tips to be used with a pipettorMicropipettors for laboratory use

#### Minipet Pipettors

| ltem  | Product Description     | UOM  |
|-------|-------------------------|------|
| 49011 | Minipet Pipettor, 5 ul  | Each |
| 49012 | Minipet Pipettor, 10 ul | Each |
| 49013 | Minipet Pipettor, 25 ul | Each |
| 49015 | Minipet Pipettor, 50 ul | Each |
| 49017 | Minpet Pipettor, 100 ul | Each |



#### StatPet<sup>™</sup> Micropipette


| Item  | Product Description                         | UOM  |
|-------|---|------|
| 49020 | StatPet <sup>™</sup> Micropipette, 1,000 ul | Each |



#### Micropipette Tips

Single use item

| ltem  | Product Description         | UOM    |
|-------|-----------------------------|--------|
| 49182 | Micropipette Tips, 200 ul   | 96/Box |
| 49183 | Micropipette Tips, 1,250 ul | 96/Box |



#### Cover Slips

- 22 mm x 22 mm
- No. 1 thickness

| ltem  | Product Description | UOM     |
|-------|---------------------|---------|
| 49001 | Cover Slips         | 100/Box |





# Microscope Slides

- Clear glass
- Ground edges
- 3" x 1"
- 1 mm 1.2 mm thickness

| Item  | Product Description        | UOM    |
|-------|----------------------------|--------|
| 49000 | Microscope Slides, 3" x 1" | 72/Box |







#### Centrifuge Vial / Tube

| ltem  | Product Description                               | UOM     |
|-------|---|---------|
| 49101 | A/V Disposable Collection Centrifuge Vial - 15 ml | 100/Box |
| 27845 | Centrifuge Tube - 50 ml                           | 50/Bag  |



# Reproduction

#### Stallion Urinary Catheter (MGi)

- Rounded tip with four side ports
- Accepts oral syringe or foley tip syringes
- Outside diameter .285"
- Length 54"
- Individually packaged and gamma irradiated or EtO exposed
- Rigid stylet for ease in insertion
- Single use item

| -     |                           |      |
|-------|---------------------------|------|
| ltem  | Product Description       | UOM  |
| 21405 | Stallion Urinary Catheter | Each |

#### VIP Impact<sup>™</sup> Specimen Shipper (MGi)

- Cost effective method of shipping fresh cooled semen and other temperature sensitive items
- Smaller than average overall dimensions resulting in less expensive shipping options
- Contains three mini phenolic packs and two specimen carriers to accommodate syringes, centrifuge tubes or bags
- Environmentally friendly insulated material (non-polystyrene)
- Outer box is a bold green color to stand out in priority shipments

| Item  | Product Description                      | UOM  |
|-------|--|------|
| 17820 | VIP Impact <sup>™</sup> Specimen Shipper | Each |



# Are You Shipping Semen?

See our Semen Extenders on Page 8.





There are various methods and instruments used to culture the uterus of any large animal. The uterine culture swab is made of calcium alginate and is used primarily in equine, however is effective in all species. Calcium alginate is a special medical grade media which is soluble in body fluids. It minimizes the risk of reddening or irritation during the collection process that may be caused by the fiber particles in other culture swabs. The swab stick is designed to minimize head loss or stick breakage in the uterus. Guarded swabs have a secure scored cap or outer body that has a secondary rigid tube within the primary shaft which protects the inner swab and shaft. All swabs are single use items, individually packaged and gamma irradiated or EtO exposed.

#### Tiegland Uterine Culture Swab with Soft Sheath

- 30" swab with outer tube and double-guarded with a plastic sheath
- Sheath is easily penetrated after passing through the cervix
- Calcium alginate fiber head for low pathogen toxicity
- Swab length is easily reduced for transport and plating
- Single use item

| Item  | Product Description                            | UOM    | 15 |
|-------|--|--------|----|
| 17705 | Tiegland Uterine Culture Swab with Soft Sheath | 25/Bag |    |

#### Tiegland Uterine Culture Swab

- 30" swab is single-guarded by an outer tube
- Calcium alginate fiber head for low pathogen toxicity
- Swab length is easily reduced for transport and plating
- Single use item

| Item  | Product Description           | UOM    |
|-------|-------------------------------|--------|
| 64400 | Tiegland Uterine Culture Swab | 25/Bag |

#### Double-Guarded Uterine Culture Swab with Calcium Alginate Tip (MG) PATENTED

- 30" swab is double-guarded with a rigid outer shell
- Calcium alginate fiber head for low pathogen toxicity
- Swab length is easily reduced for transport and plating
- Single use item

| ltem  | Product Description   | UOM    |
|-------|---|--------|
| 64300 | Double-Guarded Uterine Culture Swab with Calcium Alginate Tip | 25/Bag |

#### Bullet Tip Capped Culture Swab with Calcium Alginate Tip

- 30" swab is single-guarded by the capped outer tube
- Calcium alginate fiber head for low pathogen toxicity
- Swab length is easily reduced for transport and plating
- Once swab is extended, the hinged cap may also be used for sample collection
- Single use item

| ltem  | Product Description                                      | UOM    |
|-------|--|--------|
| 17715 | Bullet Tip Capped Culture Swab with Calcium Alginate Tip | 25/Bag |









#### Double-Guarded Cytology Brush & Culture Swab Combination (MGi)

- 31" stick double-guarded with a rigid outer shell
- Nylon brush removes easily for cytology
- Obtain cytology and/or culture samples in one easy application
- Calcium alginate fiber head for low pathogen toxicity
- Less trauma to the uterus and decreases pathogen introduction
- Easily reduced for transport and plating
- Single use item

| Item  | Product Description                                    | UOM    |
|-------|--|--------|
| 17726 | Double-Guarded Cytology Brush                          | 25/Bag |
| 17725 | Double-Guarded Cytology Brush Culture Swab Combination | 25/Bag |

#### Uterine Speculum (MGi)

Speculums are used to protect culture swabs, pipettes and any instrument used in uterine procedures, preventing vaginal contamination during entry. Each speculum has an extra heavy wall to prevent collapse with a rounded end to prevent insertion trauma. Individually packaged and available non-sterile or gamma irradiated.

- 18" long x 1-3/8" outer diameter
- Silver inside and out
- Rounded end is lightly lubricated
- Single use item

| Item  | Product Description        | UOM    |
|-------|----------------------------|--------|
| 63510 | Speculum, Gamma Irradiated | 50/Box |

# - EMCARE<sup>™</sup> Embryo Transfer Products

For successful embryo transfer, embryos must be able to survive outside of the body in solutions closely mimicking normal physiological conditions. The EMCARE<sup>™</sup> product

line has been developed to safely provide optimum conditions for embryos at every stage in the process of assisted reproduction. Designed for specialist veterinarians and animal reproduction researchers in the embryo transfer industry, EMCARE<sup>™</sup> media provides ideal conditions for flushing, handling, freezing and thawing embryos, with formulas developed with the same basic composition to reduce environmental stress on the embryos. EMCARE<sup>™</sup> media is pre-mixed, ready-to-use media that

pre-mixed, ready-to-use media that saves time and guarantees consistent quality.



Cytology Brush & Uterine Speculum







#### Single Line Uterine Lavage Tubing

- Length 60"
- Large diameter or conventional sized tubing with foley adapter at one end and a non-vented spike at the opposite end to penetrate media bags or bottles
- Attaches to a large bore pipette or foley catheter
- Individually packaged and gamma irradiated or EtO exposed
- Available in two diameters, high and regular flow
- Single use item

| Item  | Product Description                                       | UOM  |
|-------|---|------|
| 17410 | Single Line Lavage Tubing, Rapid Flow300" Inside Diameter | Each |
| 63625 | Single Line Lavage Tubing201" Inside Diameter             | Each |

#### Uterine Flushing Catheter

- Length 72"
- Large bore tubing with side eyelets and open end to enable complete flushing
- Non-vented spike to connect to a bag
- Halkey-Roberts<sup>®</sup> clamp to regulate flow rate
- Individually packaged and gamma irradiated or EtO exposed
- Single use item

| Item  | Product Description       | UOM  |
|-------|---------------------------|------|
| 21410 | Uterine Flushing Catheter | Each |

#### Simplex Uterine Flushing Catheter

- Length 75"
- Large bore tubing for rapid flow
- Large bore tubing with open end to enable complete flushing
- Vented funnel for large animal fluid administration
- Funnel opening is 2.6 cm inside diameter that stretches to fit over most common size bottle necks
- Halkey-Roberts<sup>®</sup> clamp to regulate flow rate
- Individually packaged and EtO exposed
- Single use item

| Item  | Product Description               | UOM  |
|-------|-----------------------------------|------|
| 16300 | Simplex Uterine Flushing Catheter | Each |

#### Simplex Funnel with Foley Adapter

- Length 72"
- Large bore tubing for rapid flow
- Vented funnel for large animal fluid administration
- Funnel opening is 2.6 cm inside diameter that stretches to fit over most common size bottle necks
- Foley adapter
- Halkey-Roberts<sup>®</sup> clamp to regulate flow rate
- Individually packaged and gamma irradiated or EtO exposed
- Single use item

| Item  | Product Description               | UOM  |
|-------|-----------------------------------|------|
| 39020 | Simplex Funnel with Foley Adapter | Each |









HALKEY-ROBERTS is a registered trademark of Halkey-Roberts Coporation.

**Global Locations:** 





Reproduction

14



#### Silicone Foley Catheters & Accessories (Mgi

- Used for flushing / irrigation, medication administration and embryo collection / transfer
- Inflatable balloon
- Silicone resists deterioration and is completely inert
- Soft and pliable
- Individually packaged and gamma irradiated

| Item  | Product Description   | UOM  |
|-------|---|------|
| 21360 | Silicone Flushing Catheter, 28FR, 75 cc Balloon, 2 Eyelets, Bullet Tip End - 135 cm | Each |
| 21362 | Silicone Flushing Catheter, 34FR, 100 cc Balloon, 2 Eyelets, Open End - 80 cm       | Each |
| 03940 | Stepped Adapter / Tubing Connector  | Each |



#### **Y-Junctions**

- Ideal for embryo transfer and uterine flushing / lavage
- 4" inflow tube connects to either a Foley or Rusch catheter
- Constructed of medical grade PVC tubing

(MQi)

- Star design tubing reduces kinking and works well in cold weather •
- Individually packaged and gamma irradiated or EtO exposed
- Single use item

#### Short Y-Junction with Foley Adapter & Syringe Connector

| Item  | Product Description                       | UOM  |
|-------|---|------|
| 04130 | Short Y-Junction with Foley Adapter - 30" | Each |

#### Long Y-Junction with Foley Adapter & Spike Connector

| Item  | Product Description                      | UOM  |
|-------|--|------|
| 04132 | Long Y-Junction with Foley Adapter - 60" | Each |

#### **Rapid Flow Y-Junction with Foley Adapter & Spike Connector**

Large bore tubing

| <ul> <li>Large bore spike for rapid flow</li> </ul> |  |      |
|---|--|------|
| Item  | Product Description                            | UOM  |
| 24010   | Rapid Flow Y-Junction with Foley Adapter - 60" | Each |

#### Long Y-Junction with Rusch Adapter & Spike Connector

| Item  | Product Description                    | UOM  |
|-------|--|------|
| 04134 | Long Rusch Y-Junction with Spike - 60" | Each |





Reproduction















| Short Y-Junction with Rusch Adapter & Syringe Connector |   |      |
|---|---|------|
| ltem  | Product Description                               | UOM  |
| 04131   | Short Rusch Y-Junction with Syringe Adapter - 30" | Each |



| Ó |  |
|---|--|
| M |  |
| J |  |

#### Long Y-Junction with Foley Adapter & Large Bore Spike Connector

| Item  | Product Description   | UOM  |
|-------|---|------|
| 04133 | Long Y-Junction with Foley Adapter & Large Bore Spike - 60" | Each |

| Syringe Flush with Foley Connector - Sh | hort Mai |
|---|----------|
|---|----------|

| Item  | Product Description                        | UOM  |
|-------|--|------|
| 24017 | Syringe Flush with Foley Connector - Short | Each |

| Syringe | Flush with Foley Connector - Long         |      |  |
|---------|---|------|--|
| Item    | Product Description                       | UOM  |  |
| 24018   | Syringe Flush with Foley Connector - Long | Each |  |

| Extensio | on Set with Syringe Adapter & Spike | (MQI) |
|----------|-------------------------------------|-------|
|          |                                     |       |

| Item  | Product Description                             | UOM  |
|-------|---|------|
| 24050 | Extension Set with Syringe Adapter & Spike, 16" | Each |

#### Syringe to Catheter Connector

Syringe Flush with Foley

| Item  | Product Description                | UOM  |
|-------|------------------------------------|------|
| 24052 | Syringe to Catheter Connector, 16" | Each |

| 3 Way Valve |       | alve (MGi)   |
|-------------|-------|--|
|             | ltem  | Product Description                                      |
|             | 24055 | 3 Way Valve - Dual Check Valve with Female Luer Lock Top |

| t |
|---|







|   | UOM  |
|---|------|
| & | Each |
|   |      |
|   |      |

Spin-Lock







#### Flush Bag with Y-Junction

- Ideal for uterine flushing / lavage
- Zipper closure with wide opening makes filling easy
- Long Y-Junction with foley adapter connected to bag
- 4" inflow tube connects to foley catheter
- Individually packaged and gamma irradiated or EtO exposed

| Single use item |   |      |
|-----------------|---|------|
| Item            | Product Description                     | UOM  |
| 04124           | Flush Bag with Y-Junction - 2-1/2 Liter | Each |
| 24100           | Flush Bag with Y-Junction - 6 Liter     | Each |



#### Equine Artificial Insemination & Infusion Kits

These kits include a sleeve and an insemination pipette, with or without a syringe, and are used in artificial insemination of a mare. Individually packaged and gamma irradiated. **Infusion Kit** includes: one 21" insemination pipette with flexible syringe adapter and one sleeve

Artificial Insemination Kit includes: one 21" insemination pipette with flexible syringe adapter, one sleeve and one all plastic, non-spermicidal 60 cc syringe

Single use item

•

| Item  | Product Description                                | UOM    |
|-------|--|--------|
| 19524 | Infusion Kit with Yellow Sof <sup>™</sup> Sleeve   | 10/Bag |
| 19525 | Infusion Kit with Ultra-Sens <sup>™</sup> Sleeve   | 10/Bag |
| 16010 | Infusion Kit with Shyr <sup>™</sup> Pro-Max Sleeve | 10/Bag |
| 19526 | Artificial Insemination Kit with 60 cc Syringe     | 10/Bag |

#### Deep Chamber Artificial Insemination Catheters

- Deep horn insemination
- Flexible shaft to follow curve of the horn
- Finger-grip control at plunger end
- STIP<sup>™</sup> Tip bulb end
- Ease in passing through cervix
- Length 29-3/4"
- Individually packaged and gamma irradiated or EtO exposed
- Single use item

| ltem  | Product Description                     | UOM  |   |
|-------|---|------|---|
| 27525 | Deep Chamber A.I. Catheter, 1 Exit Port | Each |   |
| 27526 | Deep Chamber A.I. Catheter, 2 Exit Port | Each |   |
| 27527 | Deep Chamber A.I. Catheter, 4 Exit Port | Each | 1 |
| 27840 | Extraction Stylet for 1/2 cc Straws     | Each |   |

#### Equi-Flex Pipette with Tube Protector

- Used for deep horn insemination
- Flexible 30" insemination pipette
- Rigid pipette protector with die cut end cap reduces vaginal contamination
- Individually packaged and EtO exposed
- Single use item

|       | -   |        |
|-------|---|--------|
| Item  | Product Description                         | UOM    |
| 31020 | Equi-Flex Pipette with Tube Protector - 30" | 25/Bag |



Infusion Pipettes & Pipette Protectors









#### **Infusion Pipettes**

Pipettes are made of PVC, designed to deliver medicants, semen or other materials through the cervix and into the uterus of an animal.

#### Large Bore Flex Tip Infusion Pipettes

- Accepts oral syringe or foley tip syringes
- Inside diameter is twice the size of a regular pipette
- Individually packaged and gamma irradiated or EtO exposed; bulk non-sterile
- Single use item

| Ite | m   | Product Description  | UOM    |
|-----|-----|--|--------|
| 190 | 058 | Infusion Pipette, Flex Tip - 25" Large Bore, Bulk, Non-Sterile               | 25/Bag |
| 190 | 059 | Infusion Pipette, Flex Tip - 25" Large Bore, Gamma Irradiated or EtO Exposed | 25/Bag |

## Flex Tip Infusion Pipettes

The syringe adapter is a friction fit, flexible PVC tube attached to a pipette and accepts any luer

lock or luer slip syringe.

- Rounded and polished entry tip to minimize tissue trauma
- Individually packaged and gamma irradiated or EtO exposed; bulk non-sterile
- Single use item

| Item  | Product Description   | UOM    |
|-------|---|--------|
| 19031 | Infusion Pipette, Flex Tip - 18", Bulk, Non-Sterile               | 25/Bag |
| 19051 | Infusion Pipette, Flex Tip - 21", Bulk, Non-Sterile               | 25/Bag |
| 19052 | Infusion Pipette, Flex Tip - 25", Bulk, Non-Sterile               | 25/Bag |
| 19033 | Infusion Pipette, Flex Tip - 18", Gamma Irradiated or EtO Exposed | 25/Bag |
| 19053 | Infusion Pipette, Flex Tip - 21", Gamma Irradiated or EtO Exposed | 25/Bag |
| 19057 | Infusion Pipette, Flex Tip - 25", Gamma Irradiated or EtO Exposed | 25/Bag |



#### Drilled Pipettes

The syringe adapter is drilled into the PVC tube and accepts any luer slip or luer lock syringe.

- Individually packaged and gamma irradiated or EtO exposed; bulk non-sterile
- Single use item

| ltem  | Product Description  | UOM    |
|-------|--|--------|
| 19030 | Infusion Pipette, Drilled Tip - 18", Bulk, Non-Sterile               | 25/Bag |
| 19050 | Infusion Pipette, Drilled Tip - 21", Bulk, Non-Sterile               | 25/Bag |
| 19035 | Infusion Pipette, Drilled Tip - 18", Gamma Irradiated or EtO Exposed | 25/Bag |
| 19055 | Infusion Pipette, Drilled Tip - 21", Gamma Irradiated or EtO Exposed | 25/Bag |

#### Pipette Protectors (MGi)

Protector is used to protect pipette and A.I. equipment from vaginal contamination.

- Two-piece construction with cored tube at one end and scored cap at the opposite end
- Pipette and A.I. equipment penetrate through scored end
- Individually packaged and gamma irradiated
- Single use item

| ltem  | Product Description           | UOM    |
|-------|-------------------------------|--------|
| 19065 | Pipette Protector, Clear - 8" | 25/Bag |
| 19060 | Pipette Protector, Blue - 8"  | 25/Bag |
| 19064 | Pipette Protector, Blue - 10" | 25/Bag |
| 19061 | Pipette Protector, Blue - 12" | 25/Bag |
| 19063 | Pipette Protector, Blue - 15" | 25/Bag |



Reproduction









•

• •

• •

Item 23138

#### Em Con<sup>™</sup> Filter Mai

- 75 micron stainless steel screen at base of filter
- Includes drain tube for controlling level of fluid in filter
- Individually packaged and gamma irradiated

Sure Flush<sup>™</sup> Embryo Transfer Filter 74 micron nylon mesh filter

Vented single step system

**Product Description** 

Sure Flush<sup>™</sup> E.T. Filter

Single use item

Single use item

| Item  | Product Description        | UOM  |
|-------|----------------------------|------|
| 04135 | Em Con <sup>™</sup> Filter | Each |

(MQI)





# Embryo Transfer Filters

#### MAI High Volume Embryo Transfer Filter (MGi)

Includes drain tube for controlling level of fluid in filter

Individually packaged and gamma irradiated

- 64 micron raised, cylindrical nylon filter •
- Contoured base for embryos to settle in center of dish
- Upon completion of flush, reservoir will maintain approximately 120 ml of media •
- Individually packaged and gamma irradiated
- Single use item

| Item  | Product Description         | UOM  |
|-------|-----------------------------|------|
| 04136 | MAI High Volume E.T. Filter | Each |





#### MAI Low Volume Embryo Transfer Filter (MGi)

- 64 micron cylindrical nylon filter
- Contoured base for embryos to settle
- Base may be used as searching dish
- Upon completion of flush, reservoir will maintain approximately 50 ml of media
- Individually packaged and gamma irradiated
- Single use item

| Item  | Product Description        | UOM  |
|-------|----------------------------|------|
| 04137 | MAI Low Volume E.T. Filter | Each |





UOM

Each





#### IVF Oocyte Aspiration Dish

- Double well aspiration receiver
- Two 25 ml chambers
- Lid designed to rotate 360°
- Ability to have two separate and unique collections
- Ability to place warm media into void to maintain temperature
- Single use item

| ltem  | Product Description        | UOM  |
|-------|----------------------------|------|
| 26500 | IVF Oocyte Aspiration Dish | Each |



#### Embryo Transfer Filter

- Angled filter along side wall
- Gridded bottom for embryo searching
- Individually packaged and gamma irradiated
- Single use item

| ltem  | Product Description         | UOM  |
|-------|-----------------------------|------|
| 25501 | E.T. Filter                 | Each |
| 25503 | E.T. Filter with Y-Junction | Each |



#### Em Con<sup>™</sup> IVF Oocyte Filters (MGi)

Designed for the aspiration, washing and filtering of oocytes in IVF procedures

- Available for 2 mm (Blue) and 3 mm (Red) needle tubing
- Color-coded filter caps prevent filter mix-up
- Cap has 2 inlet ports to connect the filter to the vacuum pump and the needle tubing
- Seal tight compression fitting for needle tubing prevents vacuum leakage
- Disposable, one time use filter provides a cost effective way to eliminate chemical and biological contamination between donors
- Economical, consistent, high quality plastic construction
- 75 micron stainless steel screen at base of filter
- Individually packaged and gamma irradiated
- Accessories available:
  - 2 mm and 3 mm needle tubing
  - Filter warmer with viewing window
- Single use item

| Item  | Product Description  | UOM  |
|-------|--|------|
| 26505 | Em Con <sup>™</sup> IVF Oocyte Filter for 3 mm Needle Tubing | Each |
| 26507 | Em Con <sup>™</sup> IVF Oocyte Filter for 2 mm Needle Tubing | Each |



Reproduction








## IVF Threaded Centrifuge Tube Cap

The IVF Threaded Centrifuge Tube Cap is designed for Ovum Pick-Up procedures using a centrifuge tube for oocyte collection.

- - Compression fitting provides tight and secure connection for needle tubing
  - Enhances biosecurity
  - Cost effective
  - Eliminates cleaning time and expense
  - Prevents toxic residue left behind from cleaning / sterilization process
  - Fits standard thread on commercial 50 ml centrifuge tubes (centrifuge tube not included) Single use item

| Item  | Product Description   | UOM  |
|-------|---|------|
| 26510 | IVF Threaded Centrifuge Tube Cap with 3 mm Needle Tubing Vacuum | Each |
|       | Port - Red  |      |
| 26512 | IVF Threaded Centrifuge Tube Cap with 2 mm Needle Tubing Vacuum | Each |
|       | Port - Blue   |      |











Embryo Transfer Supplies

## Em Con<sup>™</sup> Embryo Filter Holder (MQi)

Stainless steel ring provides solid support of embryo filter during embryo / oocyte collections and also allows the operator to maintain the filter at various heights.

- Stainless steel •
- Attaches to ring stand rods up to 7/8" in diameter •
- Fits both Em Con<sup>™</sup> Filter (#04135) and Em Con<sup>™</sup> IVF Oocyte Filters (#26505 & #26507)

| Item  | Product Description                      | UOM  |
|-------|--|------|
| 25801 | Em Con <sup>™</sup> Embryo Filter Holder | Each |

#### (Mgi Acrodisc<sup>®</sup> Syringe Filter

The Acrodisc® Syringe Filter is used to filter serum and media in preparation for specific procedures. Fits luer lock and luer slip syringes.

- Individually packaged and sterile
- Single use item

| Item  | Product Description                               | UOM    |
|-------|---|--------|
| 04701 | Acrodisc <sup>®</sup> Syringe Filter, 0.20 Micron | 50/Box |

#### Embiotic<sup>™</sup> III (MQi)

- Lyophilized newborn calf serum with Penicillin / Streptomycin used as an additive
- Contains 10 ml newborn calf serum, 10,000 units/ml Penicillin and 10,000 mcg/ml Streptomycin
- Vacuum sealed lyophilized product

| Item  | Product Description                    | UOM  |
|-------|--|------|
| 04430 | Embiotic <sup>™</sup> III - 10 ml Vial | Each |

#### **ICD**<sub>bio</sub> Albumin Solution

- Bovine Serum Albumin is of New Zealand origin, virus-screened and irradiated •
- Water is highly purified with a low endotoxin specification
- Albumin functions as a surfactant, energy source and detoxifier for embryos during embryo transfer procedures
- Sterile-filtered

| -     |  |       |
|-------|--|-------|
| ltem  | Product Description  | UOM   |
| 22502 | Albumin Solution, 20% w/v, 20 ml - Non-Returnable, Refrigerate | 5/Box |

ACRODISC is a registered trademark of Pall Corporation.



21



#### EMCARE<sup>™</sup> Holding Solution

- Cpbio
- For culturing embryos prior to transfer or freezing

For freezing embryos for direct embryo transfer

- Contains New Zealand Bovine Serum Albumin
- Optimized nutrient composition improves embryo survival rate
- Available in three sizes

| Item  | Product Description  | UOM    |
|-------|--|--------|
| 22506 | EMCARE <sup>™</sup> Holding Solution, 6 ml - Non-Returnable, Refrigerate   | 20/Box |
| 22508 | EMCARE <sup>™</sup> Holding Solution, 20 ml - Non-Returnable, Refrigerate  | 5/Box  |
| 22510 | EMCARE <sup>™</sup> Holding Solution, 100 ml - Non-Returnable, Refrigerate | 5/Box  |

**Ch**pio





UOM

20/Box

5/Box

# EMCARE<sup>™</sup> Glycerol 10% v/v

EMCARE<sup>™</sup> Ethylene Glycol 1.5M

**Product Description** 

• Glycerol is a cryoprotectant that must be removed from the frozen embryo before it is transferred into a recipient

EMCARE<sup>™</sup> Ethylene Glycol 1.5M, 6 ml - Non-Returnable, Refrigerate

EMCARE<sup>™</sup> Ethylene Glycol 1.5M, 20 ml - Non-Returnable, Refrigerate

| Item  | Product Description  | UOM    |
|-------|--|--------|
| 22520 | EMCARE <sup>™</sup> Glycerol 10% v/v, 6 ml - Non-Returnable, Refrigerate | 20/Box |

## EMCARE<sup>™</sup> Sucrose 34.2% w/v 1M

• For thawing embryos where Glycerol is the cryoprotectant

| Item  | Product Description  | UOM    |
|-------|--|--------|
| 22525 | EMCARE <sup>™</sup> Sucrose 34.2% w/v 1M, 6 ml - Non-Returnable, Refrigerate | 20/Box |

## EMCARE<sup>™</sup> Thawing System

- For thawing embryos where Glycerol is the cryoprotectant
- Each kit contains 4 tubes, 6 ml each. Each kit contains one of each solution:
  - Solution #1: EMCARE<sup>™</sup> 6% Glycerol & 10.3% w/v Sucrose
  - Solution #2: EMCARE<sup>™</sup> 3% Glycerol & 10.3% w/v Sucrose
  - Solution #3: EMCARE<sup>™</sup> 10.3% w/v Sucrose
  - Solution #4: EMCARE<sup>™</sup> Holding Solution

| Item  | Product Description  | UOM        |
|-------|--|------------|
| 22530 | EMCARE <sup>™</sup> Thawing System, 4 x 6 ml - Non-Returnable, Refrigerate | 4 Kits/Box |

# EMCARE<sup>™</sup> Complete Ultra Flushing Solution

- Contains New Zealand Bovine Serum Albumin providing surfactant and embryo protective properties with binding sites absorbing plasticizers, endotoxins and proteases
- MOPS buffered flushing solution
- Contains Kanamycin Sulfate and New Zealand Bovine Serum Albumin
- Media is sterile-filtered

**Global Locations:** 

Patented

| ltem  | Product Description  | UOM  |
|-------|--|------|
| 22060 | EMCARE <sup>™</sup> Complete Ultra Flushing Solution, 1 Liter Bag - Non-Returnable | Each |
| 22062 | EMCARE <sup>™</sup> Complete Ultra Flushing Solution, 2 Liter Bag - Non-Returnable | Each |







•

•

**Item** 22512

22514

22







#### **ICU** bio **EMCARE<sup>™</sup> BioFree Flushing Solution**

- Contains no items of biological origin •
- Contains PVA as a substitute for Bovine Serum Albumin
- PVA exhibits the same surfactant activity without the risk of introducing adventitious viral agents or other biochemical contaminants
- Contains Kanamycin Sulfate and PVA
- MOPS buffered flushing solution
- Media is sterile-filtered

| ltem  | Product Description   | UOM  |
|-------|---|------|
| 22065 | EMCARE <sup>™</sup> BioFree Flushing Solution, 1 Liter Bag - Non-Returnable | Each |
| 22067 | EMCARE <sup>™</sup> BioFree Flushing Solution, 2 Liter Bag - Non-Returnable | Each |

#### Mgi **Bag & Bottle Warmers**

- Flushing media stays warm in cold weather and cool in warmer temperatures
- Lightweight and durable
- Economical
- 1 & 2 liter bag warmers include carabiner and 2 Mini-gel Packs
- 6 liter bag warmer includes carabiner and 4 Mini-gel Packs

| ltem  | Product Description       | UOM   |
|-------|---------------------------|-------|
| 22808 | Warmer for 500 cc Bottle  | Each  |
| 22810 | Warmer for 1 Liter Bottle | Each  |
| 22807 | Warmer for 6 Liter Bag    | Each  |
| 22805 | Warmer for 2 Liter Bag    | Each  |
| 22800 | Warmer for 1 Liter Bag    | Each  |
| 26054 | Mini-gel Pack 4 oz        | 2/Bag |
| 26055 | Mini-gel Pack 4 oz        | 4/Bag |





# Flushing media and embryos stay warm in cold weather and cool in warmer temperatures

- Protects embryos from harmful UV light
- Windows for monitoring fluid depth
- Item #25803 has open top, pulled cord bottom. All other warmers have zippered tops
- Lightweight and durable

Filter Warmers

Economical

•

| Item  | Product Description  | UOM  |
|-------|--|------|
| 25803 | Warmer for Em Con <sup>™</sup> IVF Oocyte Filter or Em Con <sup>™</sup> Filter | Each |
| 25804 | Warmer for MAI Low Volume E.T. Filter  | Each |
| 25805 | Warmer for MAI High Volume E.T. Filter   | Each |
| 25810 | Warmer for E.T. Filter   | Each |
| 25811 | Warmer for Sure Flush <sup>™</sup> E.T. Filter                                 | Each |

## Embryo Transfer & Artificial Insemination Gun Warmer

- Keeps E.T. and A.I. gun warm in cold weather and cool in warmer temperatures
- Convenient shoulder strap
- Lightweight and durable
- Economical
- Includes 4 Mini-gel Packs

| Item  | Product Description    | UOM   |
|-------|------------------------|-------|
| 26050 | E.T. & A.I. Gun Warmer | Each  |
| 26055 | Mini-gel Pack 4 oz     | 4/Bag |





Reproduction





#### (MQi Square Gridded Dish

- Two-piece plastic square gridded search dish, 3"
- 100 mm x 100 mm x 15 mm
- Lid included
- Sterile
- Single use item

| -     |   |        |
|-------|---|--------|
| Item  | Product Description                                   | UOM    |
| 04041 | Square Gridded Dish with Lid, 100 mm x 100 mm x 15 mm | 10/Bag |





- Conical bowl-shaped polystyrene dish
- Used to search retained filtered flush media for embryos or oocytes
- Total volume capacity is 120 ml
- Gamma irradiated
- Single use item

| ltem  | Product Description | UOM    |
|-------|---------------------|--------|
| 04010 | Ova Transfer Dish   | 12/Bag |

# Four Well Dish

- Two-piece polystyrene dish for holding, washing or freezing embryos •
- Key-wayed lid has frosted area for plate coding
- Each well is numbered
- Ideal for four-step thawing of frozen embryos
- Large base prevents tipping
- Individually packaged and gamma irradiated
- Single use item

| Item  | Product Description | UOM    |
|-------|---------------------|--------|
| 04056 | Four Well Dish      | 10/Bag |

## Double Well Dish

- Two-piece polystyrene dish for holding, washing or freezing embryos •
- . Curved side of wells has a 20° radius
- Large base prevents tipping and key-wayed lid has frosted area for plate coding
- Each well is numbered
- Individually packaged and gamma irradiated
- Single use item

| Item  | Product Description | UOM    |
|-------|---------------------|--------|
| 04055 | Double Well Dish    | 20/Bag |

#### (MQI) Twelve Well Trypsin Dish

- Designed for in vivo and in vitro wash procedures
  - Twelve 2.5 ml volume wells
- Two - 1.2 ml volume rinse / discard wells
- Vented for gas exchange, open and closed positions
- Individually packaged and gamma irradiated •
- Single use item

•

| ltem  | Product Description                     | UOM    |
|-------|---|--------|
| 25150 | Twelve Well Trypsin Dish, 96 mm x 96 mm | 10/Bag |











**Embryo Transfer Supplies** 

Reproduction









## Petri Dish

- Two-piece polystyrene search dish
- 10 mm x 35 mm
- Gamma irradiated
- Single use item

| ltem  | Product Description        | UOM    |
|-------|----------------------------|--------|
| 04020 | Petri Dish - 10 mm x 35 mm | 20/Bag |

#### ESP<sup>™</sup> Tissue Culture Dishes (MGi)

- Two-piece polystyrene culture dish
- Each well is numbered

#### • Individually packaged and gamma irradiated

- ESP<sup>™</sup> Tissue Culture Dish, Square Grid
  - Single use item

| Item  | Product Description                               | UOM    |
|-------|---|--------|
| 04042 | ESP <sup>™</sup> Tissue Culture Dish, Square Grid | 10/Bag |

#### ESP<sup>™</sup> Tissue Culture Dish, Six Well

• Six - 2.5 ml volume wells

| Single use item |  |        |
|-----------------|--|--------|
| ltem            | Product Description                            | UOM    |
| 04045           | ESP <sup>™</sup> Tissue Culture Dish, Six Well | 10/Bag |

#### **ESP<sup>™</sup> Tissue Culture Dish, Micro Six Well**

- Six 1 ml volume wells
- Single use item

| Item  | Product Description                                  | UOM    |
|-------|--|--------|
| 04046 | ESP <sup>™</sup> Tissue Culture Dish, Micro Six Well | 10/Bag |

#### ESP<sup>™</sup> Tissue Culture Dish, Four Well

| • S   | ingle use item                                  |        |
|-------|---|--------|
| Item  | Product Description                             | UOM    |
| 04044 | ESP <sup>™</sup> Tissue Culture Dish, Four Well | 10/Bag |

#### ESP<sup>™</sup> Tissue Culture Dish, Well in Well

Single use item

| Item  | Product Description                                | UOM    |
|-------|--|--------|
| 04043 | ESP <sup>™</sup> Tissue Culture Dish, Well in Well | 10/Bag |

#### In Vitro Fertilization & Transfer Catheters

Designed to manipulate and/or move embryos in a dish or for dish to dish transfer.

#### **IVF Catheters**

• Inside diameter .025"

| • | Single use item |  |
|---|-----------------|--|
|---|-----------------|--|

| Item  | Product Description            | UOM  |
|-------|--------------------------------|------|
| 04183 | IVF Catheter, Non-Sterile      | Each |
| 04185 | IVF Catheter, Gamma Irradiated | Each |
|       | IVF Catheter, Gamma Irradiated | Eac  |

#### **Tom Cat Catheter**

- Inside diameter .018"
- Single use item

| ltem  | Product Description                       | UOM  |
|-------|---|------|
| 04186 | Tom Cat Catheter, 3.5FR - 4-1/2", Sterile | Each |

#### Embryo Transfer Catheter with Flex Adapter

- Fits luer slip tip syringe
- Inside diameter .022"
- Single use item

| -     |  |      |
|-------|--|------|
| Item  | Product Description  | UOM  |
| 04187 | Embryo Transfer Catheter with Flex Adapter, Gamma Irradiated | Each |



Embryo Transfer Supplies













#### Straws (MQi)

- Two straw sizes: 1/4 cc (inside diameter .060") and 1/2 cc (inside diameter .084")
- Uses a wick-and-powder method to seal contents
- Available gamma irradiated or non-sterile and in various package sizes
- All uncut straws are 5-1/4" long, all cut straws are 4-3/4" long
- Single use item

| Item   | Product Description                                     | UOM               |
|--------|---|-------------------|
| 04156  | Straw, Clear - 1/4 cc                                   | 100/Bag           |
| 04170  | Straw, Clear - 1/4 cc Bulk                              | 2,000/Bag         |
| 04158  | Straw, Clear - 1/2 cc                                   | 100/Bag           |
| 04171  | Straw, Clear - 1/2 cc Bulk                              | 2,000/Bag         |
| 04174  | Straw, Yellow - 1/4 cc Cut                              | 100/Bag           |
| 04170Y | Straw, Yellow - 1/4 cc Bulk                             | 2,000/Bag         |
| 04155  | Straw, Yellow - 1/2 cc                                  | 100/Bag           |
| 04171Y | Straw, Yellow - 1/2 cc Bulk                             | 2,000/Bag         |
| 04163  | Straw, Clear, Cut - 1/4 cc, Gamma Irradiated            | 10 Bags of 5 Each |
| 04164  | Straw, Transparent Yellow, D.T 1/4 cc, Gamma Irradiated | 10 Bags of 5 Each |
| 04166  | Straw, Clear - 1/4 cc, Gamma Irradiated                 | 10 Bags of 5 Each |
| 04168  | Straw, Clear - 1/2 cc, Gamma Irradiated                 | 10 Bags of 5 Each |

#### Adapter Plugs (MQi)

Two straw plug method, one straw contains the embryo and one straw identifies the contents.

- One end fits 1/4 cc straw, opposite end fits 1/2 cc straw
- Single use item

| Item    | Product Description                     | UOM     |
|---------|---|---------|
| 27825WH | Adapter Plug, White - 1/4 cc to 1/2 cc  | 100/Bag |
| 27825YL | Adapter Plug, Yellow - 1/4 cc to 1/2 cc | 100/Bag |



#### Sealing Rods (MQi)

One straw plug method that seals a 1/4 cc straw and allows for straw identification.

Single use item

| Item   | Product Description          | UOM     |
|--------|------------------------------|---------|
| 04154  | Sealing Rod, White - 1/4 cc  | 100/Bag |
| 04154Y | Sealing Rod, Yellow - 1/4 cc | 100/Bag |

#### Straw Adapters (MQi)

Used to connect a 1/2 cc or 1/4 cc straw to a syringe.

Single use item

| Item  | Product Description     | UOM     |
|-------|-------------------------|---------|
| 27810 | Straw Adapter - 1/4 cc  | 100/Bag |
| 27812 | Straw Adapter - 1/2 cc  | 100/Bag |
| 27814 | Universal Straw Adapter | 10/Bag  |

## Straw Cutter

Ideal for quickly and safely cutting the tips from sealed semen straws.

| Item  | Product Description | UOM  |
|-------|---------------------|------|
| 04160 | Straw Cutter        | Each |













#### CITO Thaw Units (MGi

The CITO Thaw Unit is an electronic unit to heat and maintain an accurate water temperature. This is a powered, glass-lined thermos with adjustable heat settings (factory preset to 96°F). The unit includes an accurate thaw monitor temperature indicator, fitted retrieval basket and power cord.

#### CITO Thaw Unit, 12-120-Volt AC

| Item  | Product Description                     | UOM  |
|-------|---|------|
| 18483 | CITO Thaw Unit, 12-120-Volt AC, Digital | Each |

#### CITO Thaw Unit, 12-220-Volt AC

| Item  | Product Description                     | UOM  |
|-------|---|------|
| 18485 | CITO Thaw Unit, 12-220-Volt AC, Digital | Each |



#### (MQI) Plastic Straw Tweezers

Used for removing straws from racks, goblets or thaw units.

| Item  | Product Description    | UOM  |
|-------|------------------------|------|
| 04161 | Plastic Straw Tweezers | Each |

#### Goblets (MGi)

- Plastic cylindrical holders for storing straws containing semen or embryos inside a liquid • nitrogen tank
- 10 mm size holds 12 1/4 cc or 5 1/2 cc straws
- 13 mm size holds 25 1/4 cc or 10 1/2 cc straws

| Item  | Product Description    | UOM     |
|-------|------------------------|---------|
| 04910 | Goblet, Clear - 10 mm  | 100/Bag |
| 04909 | Goblet, Yellow - 10 mm | 100/Bag |
| 04911 | Goblet, Clear - 13 mm  | 100/Bag |
| 04912 | Goblet, Yellow - 13 mm | 100/Bag |



#### Canes for Goblets (MGi)

Aluminum shafts hold two goblets in place inside a liquid nitrogen tank ٠

| Item  | Product Description   | UOM     |
|-------|-----------------------|---------|
| 04915 | Cane for 10 mm Goblet | 100/Bag |
| 04916 | Cane for 13 mm Goblet | 100/Bag |



#### Tabs for Canes (MQi)

- Identification tabs to mark or code canes
- For 10 mm & 13 mm canes

| Item    | Product Description    | UOM     |
|---------|------------------------|---------|
| 27820WH | Tabs for Canes, White  | 100/Bag |
| 27820YL | Tabs for Canes, Yellow | 100/Bag |
| 27820BL | Tabs for Canes, Blue   | 100/Bag |
| 27820GN | Tabs for Canes, Green  | 100/Bag |





Reproduction



27



#### Embryo Transfer Mini-Guns

- Smallest diameter at delivery tip is .060"
- Uses IMV Blue Metallic Tip Sheath with side delivery
- Uses 1/4 cc straw
- No need to cut straws when using Deep Chamber Mini-Gun

| Item  | Product Description   | UOM  |
|-------|---|------|
| 04805 | Embryo Transfer Mini-Gun 21", Deep Chamber, 1/4 cc for Uncut Straws | Each |
| 04811 | Embryo Transfer Mini-Gun 21", 1/4 cc for Cut or Uncut Straws        | Each |
| 04806 | Replacement Donut for Mini-Gun                                      | Each |

#### Insemination Cuozzo Guns

- Straw is retained in the gun with a sheath that envelopes the outer diameter of the gun, securing the straw in place
- Stainless steel
- Use with standard 1/4 cc or 1/2 cc straws

| Item  | Product Description   | UOM  |
|-------|---|------|
| 26000 | ESTROTECT <sup>™</sup> All-2-Mate <sup>™</sup> Universal Gun, 18" for use with 1/4 cc or  | Each |
|       | 1/2 cc Straws (Use Universal Unsplit Sheaths, Blue Adapter #26300)  |      |
| 04145 | Cuozzo Gun, Standard 1/4 cc - 18" (Use Universal Split Sheaths, Green<br>Adapter and Donut or 17-1/2" STIP <sup>™</sup> Tip Sheaths #27425)       | Each |
| 04146 | Cuozzo Gun, Standard 1/2 cc - 18" (Use Universal Split Sheaths, Green Adapter and Donut or 17-1/2" STIP <sup>™</sup> Tip Sheaths #27425)          | Each |
| 04147 | Cuozzo Gun, 1/2 cc Spiral - 18" (Use Universal Unsplit Sheaths, Blue Adapter #26300)  | Each |
| 04148 | Cuozzo Gun, 1/4 cc Spiral - 18" (Use Universal Unsplit Sheaths, Blue Adapter #26300)  | Each |
| 04150 | Cuozzo Gun, Standard 1/4 cc - 21" (Use Split Sheath and Donut, 21"<br>Side Metallic Sheath #04169 or 20-1/2" STIP <sup>™</sup> Tip Sheath #27445) | Each |
| 04152 | Dual Mate <sup>™</sup> Universal Gun, 21" for use with 1/4 cc or 1/2 cc Straws (Use 20-1/2" STIP <sup>™</sup> Tip Sheath #27445)                  | Each |
| 27802 | Replacement Donut for Cuozzo Gun  | Each |

| 40.0 |   |
|------|---|
|      |   |
| ATH  |   |
|      |   |
|      |   |
| 40-0 | - |
|      |   |
| 40-1 |   |
|      |   |
|      |   |
|      |   |
| 10   |   |

## STIP<sup>™</sup> Tip Sheaths (MQi)

- Tapered tip split sheath uses 1/4 cc or 1/2 cc straws
- Ease in passing through cervix
- Tip allows for leakproof seating with straw
- Conveniently fits most guns
- Restricted diameter exit orifice projects semen in desired direction (inside diameter at tip .028")
- Tip reduces trailing of semen during withdrawal
- Dual function, use with fresh or frozen semen
- Lightweight and flexible
- Single use item

| Item  | Product Description   | UOM    |
|-------|---|--------|
| 27425 | STIP <sup>™</sup> Tip Universal Split Sheath - 17-1/2"                      | 25/Bag |
| 27445 | STIP <sup>™</sup> Tip Universal Split Sheath - 20-1/2"                      | 25/Bag |
| 27500 | STIP <sup>™</sup> Tip Universal Split Sheath with Plastic Chemise - 20-1/2" | 25/Bag |

| Global Locations: |  |
|-------------------|--|
| Global Locations: |  |



Reproduction



#### Chemise (MGi

- Plastic sheath that covers a Cuozzo or implant gun
- Helps prevent vaginal contamination
- Single use item



| Item  | Product Description   | UOM     |
|-------|-----------------------|---------|
| 04810 | Plastic Chemise - 18" | 100/Box |
| 04808 | Plastic Chemise - 21" | 80/Box  |

## Standard, Universal & Metallic Tip Sheaths

- Outer plastic shell for a Cuozzo or implant gun retains a loaded straw within the gun
- End or side delivery for 1/4 cc or 1/2 cc straws
- Single use item

| Item  | Product Description  | UOM    |
|-------|--|--------|
| 04151 | Standard Split Sheath, 1/2 cc - 18"  | 5/Bag  |
| 04167 | Standard Split Sheath, 1/2 cc without Adapter - 18"  | 50/Bag |
| 04813 | Standard Split Sheath, 1/2 cc without Adapter - 18" Individually Packaged  | 25/Bag |
| 04165 | Universal Split Sheath with Green Adapter - 18"  | 50/Bag |
| 04159 | Universal Unsplit Sheath with Blue Adapter - 18"   | 50/Bag |
| 26300 | Universal Unsplit Sheath with Blue Adapter - 18" (for use with ESTROTECT <sup>™</sup> All-2-Mate <sup>™</sup> Universal Gun) | 50/Bag |
| 04812 | Universal Split Sheath with Green Adapter - 18" Individually Packaged  | 25/Bag |
| 04817 | Universal Unsplit Sheath with Blue Adapter - 18" Individually Packaged   | 25/Bag |
| 04169 | Metallic Tip Sheath, 1/4 cc Side Delivery, Clear - 21" (fits #04150)   | 5/Bag  |
| 04814 | Metallic Tip Sheath, 1/4 cc Side Delivery, Blue - 21" (fits #04805 & #04811 Mini-Guns)                                       | 5/Bag  |
| 27530 | Metallic Tip Sheath, 1/4 cc Side Delivery, Blue - 21" Individually Packaged, Sterile (fits #04805 & #04811 Mini-Guns)        | 5/Bag  |

#### Bovine Cervical Dilator (Mgi

The Bovine Cervical Dilator is a stainless steel instrument used in embryo transfer to dilate the

cervix in order to pass a catheter for embryo implantation collection.

- Tapered 22" length
- Stainless steel

| Item  | Product Description     | UOM  |
|-------|-------------------------|------|
| 27100 | Bovine Cervical Dilator | Each |

# "The Tickler"™

**CGS Ejaculator Kit** 

## Probes available for use on:

- **Bulls**
- Rams

See

Bucks

• Elk









#### Norm-Ject<sup>®</sup> Syringes

- All plastic syringe
- Non-spermicidal and non-embryocidal
- Eccentric (off centered) tip
- Individually packaged and sterile
- Single use item

| ltem  | Product Description                                | UOM     |
|-------|--|---------|
| 04708 | HSW Norm-Ject <sup>®</sup> Syringe, 1 cc, Sterile  | 100/Box |
| 04710 | HSW Norm-Ject <sup>®</sup> Syringe, 3 cc, Sterile  | 100/Box |
| 04712 | HSW Norm-Ject <sup>®</sup> Syringe, 5 cc, Sterile  | 100/Box |
| 04714 | HSW Norm-Ject <sup>®</sup> Syringe, 10 cc, Sterile | 100/Box |
| 04716 | HSW Norm-Ject <sup>®</sup> Syringe, 20 cc, Sterile | 100/Box |
| 04719 | HSW Norm-Ject <sup>®</sup> Syringe, 30 cc, Sterile | 50/Box  |
| 04720 | HSW Norm-Ject <sup>®</sup> Syringe, 60 cc, Sterile | 30/Box  |



#### Pro-Ject<sup>™</sup> 2-Piece Syringes - Luer Slip & Luer Lock Tips

A non-spermicidal, non-embryocidal syringe commonly used in artificial insemination (A.I.) and embryo transfer (E.T.) applications.

- Shorter stroke makes syringability easier
- Larger finger grips
- Economical for one time disposable use
- Two-piece all plastic syringe
- Individually packaged and gamma irradiated
- Available in luer slip or luer lock tips
- Single use item

|       | 6   |         |
|-------|---|---------|
| ltem  | Product Description   | UOM     |
| 60020 | Pro-Ject <sup>™</sup> 2-Piece Syringe, Luer Slip Tip, 15 cc, Gamma Irradiated | 100/Box |
| 60040 | Pro-Ject <sup>™</sup> 2-Piece Syringe, Luer Slip Tip, 35 cc, Gamma Irradiated | 50/Box  |
| 60060 | Pro-Ject <sup>™</sup> 2-Piece Syringe, Luer Slip Tip, 60 cc, Gamma Irradiated | 30/Box  |
| 60080 | Pro-Ject <sup>™</sup> 2-Piece Syringe, Luer Slip Tip, 80 cc, Gamma Irradiated | 25/Box  |
| 60021 | Pro-Ject <sup>™</sup> 2-Piece Syringe, Luer Lock Tip, 15 cc, Gamma Irradiated | 100/Box |
| 60041 | Pro-Ject <sup>™</sup> 2-Piece Syringe, Luer Lock Tip, 35 cc, Gamma Irradiated | 50/Box  |
| 60061 | Pro-Ject <sup>™</sup> 2-Piece Syringe, Luer Lock Tip, 60 cc, Gamma Irradiated | 30/Box  |
| 60081 | Pro-Ject <sup>™</sup> 2-Piece Syringe, Luer Lock Tip, 80 cc, Gamma Irradiated | 25/Box  |
|       |   |         |







# Looking for More Syringe Options?



# See Pages 163-164 & 189-190.

NORM-JECT is a registered trademark of Henke Saas, Wolf GmbH.





Reproduction

Syringes for E.T. / Semen



#### Syringe Caps

Syringe Caps attach to syringe tip to allow for safe transport and/or dispensing of medication.

| ltem  | Product Description            | UOM       |
|-------|--------------------------------|-----------|
| 02201 | Cap for Luer Lock Tip Syringe  | 1,000/Bag |
| 02202 | Cap for Luer Lock Tip Syringe  | 500/Vial  |
| 02203 | Cap for Luer Lock Tip Syringe  | 100/Vial  |
| 02204 | Cap for Luer Slip Tip Syringe  | 1,000/Bag |
| 02205 | Cap for Luer Slip Tip Syringe  | 500/Vial  |
| 02206 | Cap for Luer Slip Tip Syringe  | 100/Vial  |
| 60807 | Cap for Oral Foley Tip Syringe | 1,000/Bag |
| 60808 | Cap for Oral Foley Tip Syringe | 500/Vial  |
| 60809 | Cap for Oral Foley Tip Syringe | 100/Vial  |



1 as

0

0

0

ESC

# Syringe Accessories

# DVM Rapid Test<sup>™</sup> II

#### Multi-Test Analyzer

- Serum / Plasma IgG Bovine, Equine and Camelid
- Colostrum IgG Bovine and Equine
- Sperm Concentration Equine
- Portable
- Lightweight
- Low cost
- Proven accuracy
- Simple, one-step test procedure
- All tests may be performed in 10 minutes or less
- Direct reading of results; no charts needed

# See Page 9 for More Information.







0

85;

Θ

ESC

 $\odot$ 

256

BSC

CSC

ESE

0

 $\bigcirc$ 



## Glide Al<sup>™</sup> Lubricant (MQi

A non-spermicidal, water soluble lubricating jelly specifically formulated for use with semen

collection and reproductive procedures.

- Extra thick formulation for use with A/V liners
- Lubricant specifically formulated for semen collection
- Stays thick when warm to minimize semen and lubricant mixing together
- Non-spermicidal, non-embryocidal
- Ideal pH
- Perfect lube for collections in all species
- Water-soluble, colorless, odorless and stainless
- Irritation-free lubricant

| ltem  | Product Description                              | UOM  |
|-------|--|------|
| 17246 | Glide Al <sup>™</sup> Lubricating Jelly - 5 oz   | Each |
| 17247 | Glide AI <sup>™</sup> Lubricating Jelly - 16 oz  | Each |
| 17248 | Glide Al <sup>™</sup> Lubricating Jelly - Gallon | Each |
| 04605 | Pump for Gallon Jug                              | Each |



A non-spermicidal, water soluble lubricating jelly specifically formulated for use with obstetrical procedures, rectal exams and general use.

- Same product as the original Glide
- Slightly thinner product than Glide Al<sup>™</sup>
- Lubricant for obstetrical work, vaginal and rectal exams
- Non-spermicidal, non-embryocidal
- Ideal pH
- Water-soluble, colorless, odorless and stainless
- Irritation-free lubricant
- 5 oz tube features an easy to use flip top cap

| ltem  | Product Description                              | UOM  |
|-------|--|------|
| 17236 | Glide OB <sup>™</sup> Lubricating Jelly - 5 oz   | Each |
| 17237 | Glide OB <sup>™</sup> Lubricating Jelly - 16 oz  | Each |
| 17238 | Glide OB <sup>™</sup> Lubricating Jelly - Gallon | Each |
| 04605 | Pump for Gallon Jug                              | Each |



GLIDE

GLDE A

# CLIDE Non-spermicidal Lubricants

2 formulas to meet your needs

Global Locations:







#### Shoulder-Length Plastic Sleeves - Non-Sterile & Gamma Irradiated (MGI)

A complete selection of shoulder-length, disposable sleeves in assorted materials and thicknesses used

for various large animal procedures. All gamma irradiated sleeves are individually packaged.

#### Ultra-Sens<sup>™</sup> Sleeves

- 1 mil copolymer •
- Length 36"
- Super sensitive feel
- Clear
- Non-sterile or gamma irradiated
- Single use item

| Item  | Product Description   | UOM     |
|-------|---|---------|
| 19535 | Ultra-Sens <sup>™</sup> Sleeve - 1 mil                                    | 100/Box |
| 19532 | Ultra-Sens <sup>™</sup> Sleeve - Individually Packaged & Gamma Irradiated | 25/Bag  |

#### Yellow Sof<sup>™</sup> Sleeves

- 1 mil copolymer
- Length 36"
- Super sensitive feel
- Non-sterile or gamma irradiated
- Single use item

| Item  | Product Description   | UOM    |
|-------|---|--------|
| 32206 | Yellow Sof <sup>™</sup> Sleeve - 1 mil                                    | 50/Bag |
| 19225 | Yellow Sof <sup>™</sup> Sleeve - Individually Packaged & Gamma Irradiated | 25/Bag |

#### Shyr<sup>™</sup> Pro-Max Sleeves

- 1.25 mil polyethylene
- Length 35"
- For a large hand
- Non-sterile or gamma irradiated
- Single use item

| U     |   |         |
|-------|---|---------|
| Item  | Product Description   | UOM     |
| 19205 | Shyr <sup>™</sup> Pro-Max Sleeve, Grey - 1.25 mil                                 | 100/Bag |
| 19222 | Shyr <sup>™</sup> Pro-Max Sleeve, Grey - Individually Packaged & Gamma Irradiated | 25/Bag  |

#### Shyr<sup>™</sup> Professional Sleeves

- 1.25 mil polyethylene
- Length 35"
- For a small hand
- Non-sterile or gamma irradiated
- Single use item

| ltem  | Product Description   | UOM     |
|-------|---|---------|
| 19200 | Shyr <sup>™</sup> Professional Sleeve, Green - 1.25 mil                                 | 100/Bag |
| 19220 | Shyr <sup>™</sup> Professional Sleeve, Green - Individually Packaged & Gamma Irradiated | 25/Bag  |















#### **Poly Sens Sleeve, Blue**

- 1.25 mil polyethylene .
- Length 35"
- Non-sterile

• Single use item

| Item  | Product Description               | UOM     |
|-------|-----------------------------------|---------|
| 32216 | Poly Sens Sleeve, Blue - 1.25 mil | 100/Bag |

#### **Disposable Sleeve**

- 1.25 mil polyethylene •
- Length 35" .
- Non-sterile
- Single use item

| -     |                                     |         |
|-------|-------------------------------------|---------|
| ltem  | Product Description                 | UOM     |
| 32000 | Disposable Sleeve, Clear - 1.25 mil | 100/Box |

#### **Magnum Sleeve, Yellow**

- 1.10 mil copolymer
- Length 39" .
- Extremely sensitive feel
- Non-sterile .
- Single use item

|       | -                                |         |
|-------|----------------------------------|---------|
| Item  | Product Description              | UOM     |
| 32210 | Magnum Sleeve, Yellow - 1.10 mil | 100/Bag |

#### Magnum Sleeve, Brown

- 1.25 mil copolymer •
- Length 39" •
- High sensitivity and durability •
- Non-sterile
- Single use item

| ltem  | Product Description             | UOM     |
|-------|---------------------------------|---------|
| 32213 | Magnum Sleeve, Brown - 1.25 mil | 100/Bag |









#### Gripper<sup>™</sup> Shoulder-Length Plastic Sleeves - Non-Sterile & Gamma Irradiated

Disposable sleeve with a comfortable elastic band that retains the disposable sleeve onto the user's arm

- without the use of shrouds, hemostats or other devices. Simple yet effective for maintaining sleeve positioning.
  - Ideal for prolonged use
  - Available non-sterile, bulk or gamma irradiated, individually packaged

#### Gripper<sup>™</sup> Ultra-Sens<sup>™</sup> Sleeves

- 1 mil copolymer
- Length 36"
- Super sensitive feel
- Clear
- Non-sterile or gamma irradiated .
- Single use item

| ltem  | Product Description  | UOM    |
|-------|--|--------|
| 32100 | Gripper <sup>™</sup> Ultra-Sens <sup>™</sup> Sleeve - 1 mil                                    | 25/Bag |
| 32110 | Gripper <sup>™</sup> Ultra-Sens <sup>™</sup> Sleeve - Individually Packaged & Gamma Irradiated | 25/Bag |

Individually packaged & gamma irradiated.

Mgi

Gripper<sup>™</sup> Yellow Sof<sup>™</sup> Sleeve

- 1 mil copolymer
- Length 36"
- Super sensitive feel
- Non-sterile
- Bulk
- Single use item

| ltem  | Product Description   | UOM    |
|-------|---|--------|
| 19300 | Gripper <sup>™</sup> Yellow Sof <sup>™</sup> Sleeve - 1 mil | 25/Bag |



Shoulder-length sleeves are made of latex material for arm and shoulder protection with a latex glove glued directly to the latex sleeve. A retention strap is affixed to the sleeve's shoulder to secure the sleeve onto the user's body.

| Item  | Product Description   | UOM  |
|-------|---|------|
| 32052 | Latex OB Sleeve with Shoulder Strap - Made to order                   | Each |
| 32054 | Latex OB Sleeve & Glove with Shoulder Strap - Large - Made to order   | Each |
| 32055 | Latex OB Sleeve & Glove with Shoulder Strap - X-Large - Made to order | Each |



#### (MQI) Ultra Sens<sup>™</sup> Ultrasound Probe Cover

| Item  | Product Description                            | UOM    |
|-------|--|--------|
| 52600 | Ultra Sens <sup>™</sup> Ultrasound Probe Cover | 25/Bag |



Australia: +61 (0) 2 6372 6622





Sleeves & Ultrasound Probe Covers



#### ProBang<sup>™</sup> Uterine Baton (MQi)

#### Simplifies retracting and replacing an everted uterus in dairy and beef cows

The ProBang<sup>™</sup> Uterine Baton is an essential instrument used to manually replace a prolapsed uterus post calving in dairy and beef cows. The unique design consists of two blunted bulbous ends and a plastic or stainless steel shaft, which acts as a third arm for the practitioner while replacing the everted uterus into the patient's body cavity.

- Bulb ends are uniquely designed for size and surface area when in contact with the uterus and user
- Either bulb may be used for the uterus while the other is placed against the abdomen of the user
- Baton may be used on a patient standing or in a down position
- Plastic bulbs are welded to shaft and are solid and impervious for ease in maintenance
- Two batons available:
  - Lightweight: Plastic shaft, .85 lbs ٠
  - Heavyweight: Stainless steel shaft, 3.30 lbs

| ltem    | Product Description                              | UOM  |  |
|---------|--|------|--|
| 19002HW | ProBang <sup>™</sup> Uterine Baton - Heavyweight | Each |  |
| 19002LW | ProBang <sup>™</sup> Uterine Baton - Lightweight | Each |  |

#### Breech Retractor<sup>™</sup> (MGi)

Stainless steel instrument used to aid in the reposition of a breech presentation of a fetus.

- Length 35-1/2"
  - Blunted, rounded tip with a circular stop approximately six inches from the end is used by inserting the instrument tip into the fetus' rectum
- Circular stop on the instrument prevents over insertion and potential trauma to rectal walls. With physical force pushing the fetus anteriorly, the practitioner attains a physical advantage for repositioning the fetus for birth

| ltem  | Product Description           | UOM  |
|-------|-------------------------------|------|
| 17040 | Breech Retractor <sup>™</sup> | Each |

#### Ellis Uterine Detorsion Rod (MGi)

Stainless steel instrument used to correct torsion or rotation of the uterus that may restrict the fetus from passing though the birth canal.

- Length 36-3/4"
- Rod shaft will not twist during use
- Single eyelet end for fetal feet manipulation
- Opposite end features groove plates for attachment of OB chains

| Item  | Product Description         | UOM  |
|-------|-----------------------------|------|
| 17010 | Ellis Uterine Detorsion Rod | Each |

#### Fetal Repeller (MGi)

Stainless steel instrument used to aid in the correction of wry-neck or head-back birth position of fetus.


- Length 35-1/2" with a handle at one end and U-shaped yoke at the other end
- U-shaped yoke is placed on the shoulder opposite of the neck and head direction
- With even force placed on the instrument, a practitioner has a physical advantage in repositioning the fetus' head for delivery

| Item  | Product Description | UOM  |
|-------|---------------------|------|
| 17042 | Fetal Repeller      | Each |





Reproduction



36










#### ESTROTECT<sup>™</sup> Breeding Indicator with Breeding Bullseye<sup>™</sup> ESTROTECT<sup>™</sup> PATENTED

The newly redesigned ESTROTECT<sup>™</sup> Breeding Indicator features patented Breeding Bullseye<sup>™</sup> technology that indicates when cows are ready to breed, helps confirm pregnancies and provides pinpoint timing for A.I. protocols.

The precision-designed Breeding Bullseye<sup>™</sup> takes all of the guesswork out of your breeding program. Simply, when the Breeding Bullseye<sup>™</sup> - or the equivalent amount - is rubbed off, it is time to breed.

| Item  | Product Description   | UOM  |
|-------|---|------|
| 78010 | ESTROTECT <sup>™</sup> Breeding Indicator with Breeding Bullseye <sup>™</sup> - | Each |
|       | Red / Orange, 50 Count  |      |
| 78019 | ESTROTECT <sup>™</sup> Breeding Indicator with Breeding Bullseye <sup>™</sup> - | Each |
|       | Blue, 50 Count  |      |
| 78016 | ESTROTECT <sup>™</sup> Breeding Indicator with Breeding Bullseye <sup>™</sup> - | Each |
|       | Green, 50 Count   |      |
| 78013 | ESTROTECT <sup>™</sup> Breeding Indicator with Breeding Bullseye <sup>™</sup> - | Each |
|       | Fuchsia, 50 Count   |      |
| 78022 | ESTROTECT <sup>™</sup> Breeding Indicator with Breeding Bullseye <sup>™</sup> - | Each |
|       | Yellow, 50 Count  |      |
| 78400 | ESTROTECT <sup>™</sup> Rubber Brush   | Each |
| 78500 | ESTROTECT <sup>™</sup> Cleaning Cloth   | Each |



# Heat Detection

Reproduction

#### **ESTROTECT<sup>™</sup>** Heat Detectors **ESTROTECT<sup>™</sup>**

The ESTROTECT<sup>™</sup> Heat Detector is an economical, simple to use and highly effective heat detection aid. It features simple "scratch-off" technology for superior results in heat detection programs. Available in three fluorescent colors - red / orange, green and fuchsia. ESTROTECT<sup>™</sup> Heat Detectors are a management tool for artificial insemination breeding programs, recipients and donors, as well as natural breeding programs.

- Unequaled accuracy
- No messy glue, easy to use
- Cost effective
- Detects standing heat as good as a bull
- Aids in determining estrus cycling, puberty status and pregnancy
- Visit www.estrotect.com to view recent university studies validating the effectiveness of ESTROTECT<sup>™</sup>

| Item  | Product Description                           | UOM  |  |  |
|-------|---|------|--|--|
| 78000 | ESTROTECT <sup>™</sup> Red / Orange, 50 Count | Each |  |  |
| 78200 | ESTROTECT <sup>™</sup> Green, 50 Count        | Each |  |  |
| 78300 | ESTROTECT <sup>™</sup> Fuchsia, 50 Count      | Each |  |  |
| 78400 | ESTROTECT <sup>™</sup> Rubber Brush           | Each |  |  |
| 78500 | ESTROTECT <sup>™</sup> Cleaning Cloth         | Each |  |  |



#### Just how effective are ESTROTECT<sup>™</sup> Heat Detectors?

A 2005 South Dakota State University study showed heat detection accuracy using ESTROTECT" Heat Detector scratch-off technology to be greater than 90% comparable to vigorous visual observation (every 3 hours) and penile deviated bulls. There is simply no heat detection product that effectively shows evidence of multiple mountings other than ESTROTECT" Heat Detectors.









# Glide Al<sup>™</sup> & Glide OB<sup>™</sup> Lubricants



ltem #17246

#### GLIDE AI<sup>™</sup> LUBRICANT

Specifically formulated for use in Semen Collection

#### and Reproductive Procedures

- Thicker consistency formulated to prevent contact with collected semen in collection sleeve
- Perfect lube for collection in all species
- Non-spermicidal, non-embryocidal
- Ideal pH
- Water-soluble, colorless, odorless and stainless
- Irritation-free lubricant



Item #17236

#### GLIDE OB<sup>™</sup> LUBRICANT

Specifically formulated for General Use and

#### **Obstetrical Applications**

- Non-spermicidal, non-embryocidal
- Ideal pH
- Water-soluble, colorless, odorless and stainless
- Irritation-free lubricant

# **BOTH PRODUCTS AVAILABLE** IN **3 SIZES!**

| Item  | Product Description                                | UOM  | ltem  | Product Description                                | UOM  |
|-------|--|------|-------|--|------|
| 17246 | Glide Al <sup>™</sup> Lubricating Jelly - 5 oz     | Each | 17236 | Glide OB <sup>™</sup> Lubricating Jelly - 5 oz     | Each |
| 17247 | Glide AI <sup>™</sup> Lubricating Jelly - 16 oz    | Each | 17237 | Glide OB <sup>™</sup> Lubricating Jelly - 16 oz    | Each |
| 17248 | Glide AI <sup>™</sup> Lubricating Jelly - 1 Gallon | Each | 17238 | Glide OB <sup>™</sup> Lubricating Jelly - 1 Gallon | Each |
| 04605 | Pump for Gallon Jug                                | Each | 04605 | Pump for Gallon Jug                                | Each |

MAI Animal Health - UNITED STATES 1.800.447.0687 sales@maianimalhealth.com MAI Animal Health - EUROPE +44 (0) 1225 791 184 sales.uk@maianimalhealth.com MAI Animal Health - AUSTRALIA +61 (0) 2 6372 6622 sales.au@maianimalhealth.com MAI Animal Health - NEW ZEALAND +64 (0) 9 444 2115 sales.nz@maianimalhealth.com

