

PowerFloat Wolf Tooth Elevator Set:

Handle and extension with choice of 4 sizes of non-slip, cutting/elevators makes wolf tooth extractions easy. The 1/3 concave design of each cutting surface prevents accidental slippage from the extraction interface. After local anaesthesia, the elevators are secured on the caudal aspect of the wolf tooth. A hammer is used to gently percuss the cutting edge up into the ligamentous attachment. The long offset handle gives adequate control and leverage to elevate the wolf tooth.



Incisor extraction Kit



Gingival and root elevator kit

